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1. McArdle N, Devereux G, Heidarnejad H et al. Long-term use of CPAP therapy for sleep apnea/hypopnea syndrome. *Am J Respir Crit Care Med* 1999; 159 (4 Pt 1):1108-14.
2. Popescu G, Latham M, Allgar V et al. Continuous positive airway pressure for sleep apnoea/hypopnoea syndrome: usefulness of a two-week trial to identify factors associated with long-term use. *Thorax* 2001; 56 (9):727-33.
3. Waldhorn RE, Herrick TW, Nguyen MC et al. Long-term compliance with nasal continuous positive airway pressure therapy of obstructive sleep apnea. *Chest* 1990; 97 (1):33-8.
4. Bradley TD, Phillipson EA. Pathogenesis and pathophysiology of the obstructive sleep apnea syndrome. *Med Clin North Am* 1985; 69 (6):1169-85.
5. Deegan PC, McNicholas WT. Pathophysiology of obstructive sleep apnoea[see comments]. *Eur Respir J* 1995; 8 (7):1161-78.
6. Kribbs NB, Pack AI, Kline LR et al. Objective measurement of patterns of nasal CPAP use by patients with obstructive sleep apnea [see comments]. *Am Rev Respir Dis* 1993; 147 (4):887-95.
7. Massie CA, Hart RW, Peralez K et al. Effects of humidification on nasal symptoms and compliance in sleep apnea patients using continuous positive airway pressure [see comments]. *Chest* 1999; 116 (2):403-8.
8. Ball EM, Banks MB. Determinants of compliance with nasal continuous positive airway pressure treatment applied in a community setting. *Sleep Medicine* 2001; 2 (3):195-205.
9. Hui DS, Choy DK, Li TS et al. Determinants of continuous positive airway pressure compliance in a group of Chinese patients with obstructive sleep apnea. *Chest* 2001; 120 (1):170-6.
10. Rauscher H, Formanek D, Popp W et al. Self-reported vs measured compliance with nasal CPAP for obstructive sleep apnea. *Chest* 1993; 103 (6):1675-80.
11. Cole P. The respiratory role of the upper airways. Toronto: Mosby-Year Book Inc; 1992. 164 p.
12. Togias AG, Naclerio RM, Proud D et al. Nasal challenge with cold, dry air results in release of inflammatory mediators. Possible mast cell involvement. *J Clin Invest* 1985; 76 (4):1375-81.
13. Constantinidis J, Knobber D, Steinhart H et al. Finestructural investigations of the effect of nCPAP-mask application on the nasal mucosa. *Acta Otolaryngol* 2000; 120 (3):432-7.
14. Nino-Murcia G, McCann CC, Bliwise DL et al. Compliance and side effects in sleep apnea patients treated with nasal continuous positive airway pressure. *West J Med* 1989; 150(2):165-9.
15. Pepin JL, Leger P, Veale D et al. Side effects of nasal continuous positive airway pressure in sleep apnea syndrome. Study of 193 patients in two French sleep centers [see comments]. *Chest* 1995; 107 (2):375-81.
16. Williams R, Rankin N, Smith T et al. Relationship between the humidity and temperature of inspired gas and the function of the airway mucosa. *Crit Care Med* 1996; 24 (11):1920-9. Fisher & Paykel Healthcare **15**
17. Martins de Araujo MT, Vieira SB, Vasquez EC et al. Heated humidification or face mask to prevent upper airway dryness during continuous positive airway pressure therapy [see comments]. *Chest* 2000; 117 (1):142-7.
18. Hayes MJ, McGregor FB, Roberts DN et al. Continuous nasal positive airway pressure with a mouth leak: effect on nasal mucosal blood flow and nasal geometry. *Thorax* 1995; 50 (11):1179-82.
19. Richards GN, Cistulli PA, Ungar RG et al. Mouth leak with nasal continuous positive airway pressure increases nasal airway resistance. *Am J Respir Crit Care Med* 1996; 154 (1):182-6.
20. Pepin JL, Krieger J, Rodenstein D et al. Effective compliance during the first three months of continuous positive airway pressure. A European prospective study of 121 patients. *Am J Respir Crit Care Med* 1999; 160 (4):1124-9.
21. Weaver TE, Kribbs NB, Pack AI et al. Night-to-night variability in CPAP use over the first three months of treatment. *Sleep* 1997; 20 (4):278-83.
22. Rosenthal L, Gerhardstein R, Lumley A et al. CPAP therapy in patients with mild OSA: implementation and treatment outcome. *Sleep Medicine* 2000; 1 (3):215-20.
23. Kline LR, Carlson P. Acceptance and compliance with continuous positive airway pressure is altered by humidification (abstract). *Sleep* 1999; 22 (Supplement):S230.
24. Rakotonanahary D, Pelletier-Fleury N, Gagnadoux F et al. Predictive Factors for the Need for Additional Humidification During Nasal Continuous Positive Airway Pressure Therapy. *Chest* 2001; 119 (2):460-5.
25. Wiest GH, Fuchs FS, Brueckl WM et al. In vivo efficacy of heated and non-heated humidifiers during nasal continuous positive airway pressure (nCPAP)-therapy for obstructive sleep apnoea [In Process Citation]. *Respir Med* 2000; 94 (4):364-8.
26. Orec R, Richards G, Cornere B et al. Movement of bacterial contaminants is reduced in heated ventilator circuits. *Respir Care* 1995; 40 (11):1201.
27. Bacon JP, Farney RJ, Jensen RL et al. Nasal continuous positive airway pressure devices do not maintain the set pressure dynamically when tested under simulated clinical conditions [In Process Citation]. *Chest* 2000; 118 (5):1441-9